VERMEIREN

MAT-X1, MAT-X2, MAT-X3

MAT-X1

INSTRUCTION MANUAL
MANUEL D'UTILISATION
GEBRUIKSAANWIJZING
GEBRAUCHSANWEISUNG
MANUALE DI ISTRUZIONI
MANUAL DE INSTRUCCIONES
INSTRUKCJA OBSŁUGI
NÁVOD K OBSLUZE

MAT-X2

MAT-X3

EN

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EN

Contents

1. Technical details	2
2. Components	2
3. Location identification plate	2
4. Explanation of symbols	2
5. Check after delivery	2
6. Intended use	3
7. Grounding	3
8. Installation	3
9. Compressor panel	3
10. Operation	3
11. Safety instructions	3
12. EMC	4
13. Troubleshooting	6
14. Care	6
15. Disinfection	6
16. Repairs / Service / Inspection	6
17. Shipping and storage	7
18. Warranty	7
19. Disposal	7
20. Declaration of conformity	7

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EN INSTRUCTION MANUAL
Alternating pressure mattresses
MAT-X1, MAT-X2, MAT-X3

1. Technical details

⚠ WARNING: Risk of unsafe settings
Use only the settings described in this manual.

Make	Vermeiren
Address	Vermeirenplein 1/15, B-2920 Kalmthout
Type	 Alternating pressure mattresses
Model	MAT-X1, MAT-X2, MAT-X3
Environment	
Operate temperature	+10°C to +40°C
Operate humidity	30 - 75% RH
Storage temperature	-10°C to +60°C
Storage humidity	15 - 90% RH
Compressor	
Width	110 mm
Length	130 mm
Height	100 mm
Length of power cable	2m
Weight	1,2 kg
Power supply	AC 220V/50Hz 8W
Cycle time	5 min / 10 min
Fuse	F1,5A, 250 V_{AC}
Protection class	Class II, Type B
Inflation pressure	6 LPM
Noise level	Extremely silent
Mattress	
Width	900 mm
Length	2000 mm
Height	MAT-X1: 65 mm MAT-X2: 90 mm MAT-X3: 114 mm
Weight	MAT-X1: 2 kg MAT-X2: 2 kg MAT-X3: 6,5 kg
Maximum user weight	MAT-X1: 125 kg MAT-X2: 150 kg MAT-X3: 130 kg
Material (cover)	MAT-X1, MAT-X2: PVC MAT-X3: Nylon-TPU
We reserve the right to introduce technical changes. Measurement tolerance +/- 15 mm / 1,5 kg	

First of all we wish to thank you for the trust you placed in us by selecting a **VERMEIREN** product. Before using the product, read the instruction manual carefully: it will familiarise you with this product.

Please remember that observance of our hints would keep your product in an excellent condition of preservation and it will function perfectly even after years of use.

If you have any further questions, please consult your specialist dealer.

To find a service facility or specialist dealer near you, contact the nearest Vermeiren facility. A list of Vermeiren facilities can be found on the last page.

2. Components

3. Location identification plate

The location of the identification plate is on the retail box.

4. Explanation of symbols

Maximum mass

Indoor use

CE conformity

Observe the safety instructions

Separate recovery and recycling of electric and electronic devices

Protection class II

Type B applied part

Keep dry

5. Check after delivery

Unpack your product and verify if the delivery is complete. Following items should be included:

- mattress, air hoses
- compressor
- manual

Verify your product for transport damage. If you find any damages after delivery proceed as follows:

- contact the transporter
- have a list made of any problems
- contact your supplier

6. Intended use

The alternating pressure mattresses are placed over the standard mattress. They can be used at home, home care or in an hospital. Preferably used on hospital or nursing beds, but you can also use them on a standard bed.

It has been designed for prevention of bedsores, affordable solution to 24 hours pressure area care. The mattress is intended for prevention and treatment of pressure ulcers / decubitus, while optimizing the patient comfort.

Do not use the mattress on patient conditions for which the application of pressure relieving therapy on alternation system is contraindicated (cervical skeletal traction, unstable spinal cord injuries).

7. Grounding


Before any connection to the output connectors is made, the compressor shall be connected to a protective earth conductor with the three-core main cable; the mains plug shall be inserted not be negated by the use of an extension cable without protective conductor.

8. Installation

⚠ DANGER: Always unplug the product immediately after using.

1. Unpack the mattress and compressor.
2. Place the overlay mattress MAT-X1 or MAT-X2 or MAT-X3 onto the bed mattress. Check that the bed mattress and mattress MAT-X1 or MAT-X2 or MAT-X3 match in size. Please note the foot end of the overlay mattress.
3. Mount the compressor to the foot panel and check the compressor is properly attached.
4. Connect the air hoses to the mattress and the compressor. Check and ensure the air hoses are not kinked or tucked under the mattress.
5. Plug the power cable into the outlet socket (grounded AC socket 230V). Check the main voltage of the outlet socket is correct before inserting the plug of the power cable.

9. Compressor panel

10. Operation

⚠ DANGER: Always unplug the product immediately after using.

⚠ DANGER: Do not reach a mattress that has fallen into water. Unplug immediately.

NOTE: Users can adjust pressure of air mattress to a desired softness by themselves, or by healthcare professional's suggestion.

NOTE: If the pressure level is consistently low, check for any leakage (tubes or connecting hoses). If necessary, replace any damaged tubes or hoses or contact local qualified dealer for repair.

1. Switch on the compressor (power switch to "ON"). The power switch and normal pressure indicators light up.
2. The compressor starts pumping air into the mattress. It takes approx. 10 minutes to inflate.
3. Turn the pressure dial to "MAX". When the mattress is inflated for the first time. Afterwards can the user adjust the air mattress to the desired softness by turning on the pressure switch of the compressor.
4. Check that your hand can be between the patient and the inflated mattress without any difficulty to see the mattress is inflated to the correct pressure.
5. A cotton sheet can be placed over the mattress to avoid skin contact and for the patient's comfort. Take care to avoid any folds in the sheet.

11. Safety instructions

- Always unplug the mattress immediately after using.
- Do not use while bathing.
- Do not place or store the product where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.
- Do not reach for a mattress that has fallen into water. Unplug immediately!
- Protect the compressor from humidity or direct moisture.
- This product should never be left unattended when plugged in.
- Close supervision is necessary when this product is used near children or invalids.
- Use this product only for its intended use as described in this manual. Do not use hoses not recommended by the manufacturer.
- The product has been dropped or damaged, or dropped into water. Return the product to local retail or service center for repairing.
- Never operate this product if it has a damaged cable or the plug is not working properly.
- Keep the cable away from heated surfaces.

- Do not use the mattress system in the immediate vicinity of heat sources.
- Do not use the mattress system in the immediate vicinity of inflammable gases or in areas where there is a risk of explosion.
- Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where the openings may be blocked, keep the air opening free of lint, hair, and other similar particles.
- Keep the tubing free from kinks.
- Never drop or insert any object into any opening or hose.
- Connect this product to a properly grounded outlet only and check the main voltage is correct.
- Do not use outdoors or operates where aerosol (spray) products are being used or where oxygen is being administered.
- Do not use the mattress system in an uninflated state.
- Remove the plug when moving the bed.
- Only pull on the plug and not on the cable.
- Keep sharp objects away from the mattress.
- Do not use unsuitable methods to fix the overlay mattress.
- Use only genuine spare parts and expendables.
- Check the specified minimum clearance for the mattress when using side rails. Eventually use a height extension for the side rail.
- **The max. user weight of the product is MAT-X1: 125 kg, MAT-X2: 150 kg, MAT-X3: 130 kg**
- Only allow repairs by the manufacturer or an authorized dealer.
- Take note of the instructions for care and service. The manufacturer is not liable for damage caused by improper servicing/care.

- Mobile medium-range transmitting and receiving devices (external car antennas)
 - Fixed transmitting and receiving sets
 - Hands-free mobile speaking devices (fixed installation)
 - Fixed radio, TV and navigation systems
- Other home appliances
 - CD player
 - Laptop
 - Microwave oven
 - Cassette recorder
 - etc.

Devices such as electric shavers and hair dryers should have no influence. However they need to function properly and have proper wiring regardless of their influence. Please read the appropriate instruction manual from the relevant manufacturer.

To reduce the effect of electromagnetic sources of interference, please read the following warnings:

- ⚠ Do not use any portable television or radio devices in the immediate vicinity of the mattress.
- ⚠ Do not use any walkie talkies or mobile phones in the immediate vicinity of the mattress.
- ⚠ Should any unintended movements occur, remove the power plug for the mattress from the wall socket.

12. EMC

Disruption of electrical equipment from electromagnetic fields in the surroundings can not be entirely avoided. Possible consequences are:

In the presence of very strong or long-lasting disruptive fields, the electrical systems might break down completely and suffer permanent damage. Possible sources of radiation include:

- Portable receiver and transmission devices (receiver and transmitter with fixed antennae)
- Walkie-talkies, cell phones or cordless telephones
- portable TV, radio and navigation devices
- other personal broadcast devices

Guidelines and manufacturer's declaration - Electromagnetic transmissions		
The mattress is intended for operations in surroundings described below. The customer or mattress user needs to make sure that it will be operated in surroundings of this type.		
Electromagnetic radiation measurements	Agreement	Electromagnetic surroundings - Guidelines
HF - Radiation pursuant to CISPR 11 (partial)	<u>Group 1</u>	The mattress uses HF energy exclusively for its internal function. Therefore its HF radiation is very low and it is unlikely that neighboring electronic devices will be disrupted.
HF - Radiation pursuant to CISPR 11 (partial)	Class B	The mattress is suitable for use in all facilities which are solely used as living areas and which have an immediate connection to the public power supply, which also supplies the building which is being used as living quarters.
Transmission of harmonics in accordance with IEC 61000-3-2	Class A	
Transmission of Voltage fluctuations/ flickers pursuant to IEC 61000-3-3	In accordance with	

Guidelines and manufacturer's declaration - Electromagnetic resistance				Guidelines and manufacturer's declaration - Electromagnetic resistance			
The mattress is intended for operations in surroundings described below. The customer or the user of the mattress needs to make sure that it will be operated in surroundings of this type.				The mattress is intended for operations in surroundings described below. The customer or the user of the mattress needs to make sure that it will be operated in surroundings of this type.			
Resistance tests	IEC 60601 Test level	Compliance level	Electromagnetic surroundings - Guidelines	Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic surroundings - Guidelines
Discharge of static electricity (ESD) pursuant to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Flooring needs to be made of wood or concrete or covered with ceramic tiles If the floor is covered with synthetic material, then relative humidity needs to be at least 30%.	Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 Mhz to 2,5 GHz	3 Vrms = V1 3 V/m = E1	Portable and mobile RF communications equipment should be used no closer to any part of the mattress, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Quick stray electrical factors / bursts pursuant to IEC 61000-4-4	± 2 kV for power cables ± 1 kV for input and output cables	± 2 kV for power cables ± 1 kV for input and output cables	The quality of the power voltage needs to meet standards typical in a commercial or hospital environment.				Recommended separation distance
Surges pursuant to IEC 61000-4-5	± 1 kV Series mode voltage ± 2 kV Common-mode voltage	± 1 kV Series mode voltage ± 2 kV Common-mode voltage	The quality of the power voltage needs to meet standards typical in a commercial or hospital environment.				$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
Voltage drops, Short-term disruptions And Fluctuations in the supply voltage pursuant to IEC 61000-4-11	< 5 % VT (>95 % drop in VT) for ½ period 40 % VT (60 % drop in VT) for 5 periods 70 % VT (30 % drop in VT) for 25 periods < 5 % VT (>95 % drop in VT) for 5 seconds	< 5 % VT (>95 % drop in VT) for ½ period 40 % VT (60 % drop in VT) for 5 periods 70 % VT (30 % drop in VT) for 25 periods < 5 % VT (>95 % drop in VT) for 5 seconds	The quality of the power voltage needs to meet standards typical in a commercial or hospital environment. Should the user of the mattress require the continuous function even with the occurrence of disruptions to the power supply, it is recommended that the nursing bed receive undisturbed power from an emergency source or a battery.				Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol:
Magnetic field for the Supply frequency (50/60 Hz) pursuant to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields with the grid frequency should comply to those normal in the commercial or hospital environments.				
NOTE VT is the grid alternating voltage before applying the test level.				NOTE 1 At 80 MHz and 800 Mhz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
				a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed strength in the location in which the mattress is used exceeds the applicable RF compliance level above, the mattress should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the mattress.			
				b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the mattress.			
The mattress is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the mattress can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the mattress as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of the transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_1}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_1}\right]\sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

13. Troubleshooting

The following list will assist you with the troubleshooting of your mattress:

- Is there a power supply (correct voltage 230V)?
- Is the compressor switched "ON"?
- Is the plug of the power cable connected correctly?
- Are the air hoses to the mattress and compressor connected or connected correctly?
- Are there any kinks in the air hoses?
- Is the mattress damaged?

Contact your specialist dealer, if the problem stays. Never try to repair the mattress yourself. Fuse may only be replaced by qualified and authorized personnel. Remove the plug, in case of any damage.

14. Care

This is a reference procedure to clean and decontaminated the mattress, it is important to follow these procedures before using the mattress on patients.

Compressor

- Do not immerse or soak the compressor.
- Check for external damage and move the compressor to the cleaning area.
- Place the compressor on a work surface and spray or wipe the outside of the case with quaternary ammonium solution, **DO NOT** spray

any cleaning solution directly on the surface of the compressor.

- **DO NOT** use a Hydrocarbonate or Phenolic based cleaning solution as this may cause damage to the case, Allow the solution to incubate for 10 minutes or according to the package insert from the manufacturer.
- Wipe case with a clean cloth, etc.. make sure all areas are clean (top, both sides and bottom).
- Spray cloth with cleaning solution and clean front panel. **DO NOT** allow excess cleaning solution on front panel. (If solution gets inside damage will occur!) Allow surface to thoroughly dry after cleaning.
- After compressor is thoroughly dry after cleaning.
- After compressor is thoroughly cleaned and dried, proceed to plug in the compressor and test the compressor to see if it runs normally.

Mattress

- Brush off or wipe down all surfaces with soap and water before wetting with any liquid disinfectant.
- Brush or wipe down all surfaces with soap and water before applying with any liquid.
- After pre-soaking, the cover is rinsed through a regular cycle in a washer with no soap, then laundered with mild detergent.
- Repeat the process with the tubing set, spray, incubate, and then wipe clean,
- Allow it to thoroughly air dry, Once the inside is dry, turn it back wipe down the outside of the cover with disinfectant.
- Dry the mattress on SUNLESS area after cleaning.

15. Disinfection

Disinfections may only be undertaken by a skilled sanitation employee or a person trained by him/her. Consult your specialist dealer.

For sterilisation use only disinfectants suitable for treating varnished metal. Observe the instructions of the cleaning liquid used. Use only proved processes (for scrubbing disinfection) and disinfectants included in the list of the Robert Koch Institute (information can be obtained from: www.rki.de).

16. Repairs / Service / Inspection

Lifetime of the mattress is influenced by its use, storage, regular maintenance, servicing and cleaning.

You should use only genuine VERMEIREN replacement parts. Repairs and renovations may only be undertaken by trained persons. Make use of the services offered by your specialist dealer. He would gladly assist you further in regard to service and repairs.

Before you resume use of the mattress, it should be serviced according to the following criteria:

General:

- Completeness
- Check main power cable and plug if there are abrasions or excessive wear.
- Ensure mattress, compressor and air hoses are connected together correctly.
- Plug in the compressor and check the airflow from the hose connection port, the airflow should alternate between ports every cycle time.
- Check the air hoses if there are kinks or breaks, for replacement, please contact local agents or dealers.
- Check the compressor and make sure both power and power indicator is off when the power switch is turned off.

Low pressure:

- Check connectors between mattress and compressor, if they are disconnected, please reconnect them.
- Check air hoses. Ensure each single cell is not broken, Pressure switch is set to firm position, keep the hoses fully inflated and inspect for air leakage.
- Exam air leakage from cells. Ensure no leakage occurs. If any leakage occurs, please use our repair kits to seal leakage hole.

Every 6 months or for every new user, it should be serviced according to the following criteria:

- General review
- Cleaning
- Disinfection (for reuse)

Every 2 years:

- Replace the air filter
- Replace the timer motor
- Perform function test/function inspection

17. Shipping and storage

The shipping and storage of the mattress shall be according following instructions:

- Store in a dry place (between -10°C and +60°C).
- The relative humidity of the air should be between 15 - 90% RH.
- Provide sufficient covering or packaging to protect the mattress.
- The mattress must be stored without being subjected to strains. (Do not put too heavy parts on the mattress, not clamping between something, ...)
- Lay the mattress out flat and upside down.

18. Warranty

Excerpt from the "General Business Conditions".

(...)

5. The terms of the guarantee may differ from country to country. Consult your specialist dealer for the guarantee period for warranty claims.

(...)

The warranty excludes damage arising from structural changes to our products, insufficient maintenance, defective or improper handling or storage or the use of non-original parts.

Likewise, the warranty excludes parts or working parts subject to natural wear and tear.

19. Disposal

When disposing of the mattress, contact your local disposal centre or return the product to your specialist dealer who, after submitting it to a hygienic procedure, will be able to send it back to the manufacturer who will dispose of and recycle it correctly, separating it into its component materials.

The manufacturer is responsible for taking back and recycling the compressor and shall meet the requirements of European Directive 2002/96/EC on Waste Electrical and Electronic Equipment.

Packaging materials can be taken to disposal or recycling centers or to your specialist dealer.

20. Declaration of conformity

The manufacturer or his authorized representative :

VERMEIREN GROUP

Address :

Vermeirenplein 1/15
2920 Kalmthout
Belgium

declares under his sole responsibility that the CE marked devices :

Productgroup: Alternating pressure mattresses + compressor
Brand: Vermeiren
Type: MAT-X1, MAT-X2, MAT-X3

have been classified as class I, according to annex IX MDD 93/42/EEC, rule 12,

and is manufactured in full conformity with the European instructions below - including the latest modifications - and with the national law, that organizes this directions :

Medical devices directive MDD 93/42/EEC: 2007, Annex VII

and is in conformity with the relevant European harmonized standards:

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Productgroup: Alternating pressure mattresses + compressor

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Service registration form

This product (name):

was inspected (I), serviced (S), repaired (R) or disinfected (D):

By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:
By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:
By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:
By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:
By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:	By (stamp): Kind of work: I / S / R / D Date:

Belgium

Vermeiren N.V.

Vermeirenplein 1 / 15

B-2920 Kalmthout

phone: +32(0)3 620 20 20

fax: +32(0)3 666 48 94

website: www.vermeiren.be

e-mail: info@vermeiren.be

France

Vermeiren France S.A.

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59133 Phalempin

phone: +33(0)3 28 55 07 98

fax: +33(0)3 20 90 28 89

website: www.vermeiren.fr

e-mail: info@vermeiren.fr

Italy

Vermeiren Italia

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I-20020 Arese MI

phone: +39 02 99 77 07

fax: +39 02 93 58 56 17

website: www.vermeiren.it

e-mail: info@vermeiren.it

Poland

Vermeiren Polska Sp. z o.o.

ul. Łączna 1

PL-55-100 Trzebnica

phone: +48(0)71 387 42 00

fax: +48(0)71 387 05 74

website: www.vermeiren.pl

e-mail: info@vermeiren.pl

Czech Republic

Vermeiren ČR S.R.O.

Nádražní 132

702 00 Ostrava 1

phone: +420 596 133 923

fax: +420 596 121 976

website: www.vermeiren.cz

e-mail: info@vermeiren.cz

Germany

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D-40472 Düsseldorf

phone: +49(0)211 94 27 90

fax: +49(0)211 65 36 00

website: www.vermeiren.de

e-mail: info@vermeiren.de

Austria

Vermeiren Austria GmbH

Schärdinger Strasse 4

A-4061 Pasching

phone: +43(0)7229 64900

fax: +43(0)7229 64900-90

website: www.vermeiren.at

e-mail: info@vermeiren.at

Switzerland

Vermeiren Suisse S.A.

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3645 Gwatt (Thun)

phone: +41(0)33 335 14 75

fax: +41(0)33 335 14 67

website: www.vermeiren.ch

e-mail: info@vermeiren.ch

Spain / Portugal

Vermeiren Iberica, S.L.

Carratera de Cartellà, Km 0,5

Sant Gregori Parc Industrial Edifici

A 17150 Sant Gregori (Girona)

phone: +34 972 42 84 33

fax: +34 972 40 50 54

website: www.vermeiren.es

e-mail: info@vermeiren.es

Manufacturer

Vermeiren GROUP

Vermeirenplein 1 / 15

B-2920 Kalmthout

phone: +32(0)3 620 20 20

fax: +32(0)3 666 48 94

website: www.vermeiren.com

e-mail: info@vermeiren.com

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All rights reserved, including translation.