



Material			
No.	Plan	Article	Quantity
1	X002670	1547400	2
2	X005741	1914539	1
3	X005840	1915117	2
4	X005687	1914548	1
5	Y022659	20454..	1
6	M12 DIN125A	1001561	5
7	D38	1001611	2
8	M12 DIN985	1001684	3
9	M8x75 DIN931	1006110	2
10	M12-16x55 ISO7379	1914019	1
11	∅3x10	1915115	4
12	---	1992559	4
13	∅140x40	1914904	2
14	M12x1,25	1994605	2
15	M8	1998994	2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

BOM SIGMA REAR BALANCE ARM				2993722
				Plan No. E001093/E
Scale: 1:3	Drafted	Verified	Confirmed	
Date	08.06.2018	11.08.2025	11.08.2025	
Initials	SD	PM	MK	