



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2

CROSS TYPE ECLIPS X4 / 28-32-37

Plan No.

E000069/C

Scale:	1:3	Drafted	Verified	Confirmed
Date	18/12/2003	15/10/2025	15/10/2025	
Initials	SD	AZ	AS	



Bill of Material			
No.	Plan	Article	Quantity
1	P00104935 (28)	2909354	1
	Y004899 (32)	1533554	1
	Y004895 (37)	1533054	1
2	P00105011 (28)	2909454	1
	Y004900 (32)	1533454	1
	Y004565 (37)	1532954	1
3	M8x110 DIN912	1991711	1
4	X001119	1001554	2
5	X002133	1992424	4
6	Round PVC D6	H001112 (28)	0,56
	Round PVC D6	H001112 (32)	0.64
	Round PVC D6	H001112 (37)	0.72
7	X001016	1001682	1
8	X002619	1994209	2
9	M6 DIN985	1001680	2
10	M6 DIN125A	1000973	2
11	X002406 (32/37)	1994019	2
	P00116898 (28)	2992402	2
12	M6x55 ISO7380	1997879	2
13	X002136	1992414	4
14*	X002677	1994453	2
15*	Y005449	1546554	2
16	4,2x9,5 DIN7981	1992618	4/8*

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 2/2

CROSS TYPE ECLIPS X4 / 28-32-37

Plan No.

E000069/C

Scale:	1:3	Drafted	Verified	Confirmed
Date	18/12/2003	15/10/2025	15/10/2025	
Initials	SD	AZ	AS	

