

2993723


E001103 / E

SIGMA SE42 LIFT WITH SEAT TILT 45°				Plan No.
Sheet 1/2				E001103/ E
Scale: 11:50	Drafted	Verified	Confirmed	
Date	19/06/2018	11/12/2023	11/12/2023	
Initials	SD	PM	MK	

Material			
No.	Plan	Article	Quantity
1	X005836	G102700	2
2	X005837	G102800	2
3	X004780	1906368	2
4	X005838	1915122	2
5	Y023035	2051254	1
6	Y023057	2052800	2
7	M5 DIN125	1001550	4
8	M6 DIN125	1001564	4
9	M6x23	1003558	2
10	---	1004094	2
11*	---	1004109	2
12	M6x20 ISO7380	1006999	2
13	M6 DIN986	1008492	2
14	D4x14 DIN1481	1904991	2
15	SE42 / 45°	1914670	1
16	600N	1915813	2
17	M8x16 DIN912	1915411	4
18	M5 DIN986	1992804	4
19	---	1999963	2 x 0.96 m
20	---	1994748	2
21	---	1994750	2
22	M8 DIN985	1001682	4
23	M8x20 DIN7991	1903136	4

2993723

\* Not depicted

<b>SIGMA SE42 LIFT WITH SEAT TILT 45°</b> Sheet 2/2				Plan No.
				E001103/ E
Scale: 11:50	Drafted	Verified	Confirmed	
Date	19/06/2018	11/12/2023	11/12/2023	
Initials	SD	PM	MK	

E001103 / E