


Material			
No.	Plan	Article	Quantity
1	XP02085	GP15599	1
2	XP02086	GP15699	1
3	X005597	2025599	1
4	Y019481	1784204	1
5	M5x20 DIN912	1003907	2
6	M5x16 DIN912	1007198	2
7	---	1908818	2
8	D3x10	1908849	4
9	D22	1990648	2

EP00487 / A

B31 CRUTCH HOLDER TURIOS				Plan No.
				EP00487 / A
Scale: 1:4	Drafted	Verified	Confirmed	
Date	06/02/2020	06/02/2020	06/02/2020	
Initials	PM	AZ	TK	