



Material			
No.	Plan	Article	Quantity
1	Y012571	1681601	1
2	Y012557	1681700	1
3	Y012556	1676100	1
4	Y012555	1676204	1
5	M8x12-DIN916/45H	1004885	1
6	M10x19	1991432	1
7	M10-DIN934	1001672	1
8	M6x25	1006998	1
9	---	1418700	1
10	M8-DIN125	1003199	1
11	---	1005328	0,067 m <sup>2</sup>
12	X002018	H7443	1
13	M6x10	1992516	2
14	---	2991788	0,4 m
15	M6-DIN985	1001680	1
16	M6x20	1006999	1
17	---	2992176 (M)	1
		2992175 (L)	1
18	---	1991236	2
19	---	1009516	2 x 0,7 m
20	---	1009710	2 x 0,2 m
21	---	1009709	2 x 0,2 m

EP00099 / C

B66 STABILO   TYPE INOVYS			Plan No.	
			EP00099 / C	
Scale: 1:4	Drafted	Verified	Confirmed	
Date	20/12/2012	14/04/2017	14/04/2017	
Initials	BS	AZ	TK	