


| Material |                 |               |            |
|----------|-----------------|---------------|------------|
| No.      | Plan            | Article       | Quantity   |
| 1        | YP07351 (P)     | 2761965 (P)   | 1          |
|          | YP07350 (L)     | 2761865 (L)   | 1          |
| 2        | X001380         | 1009726       | 2          |
| 3        | X002437         | 1992699       | 2          |
| 3A       | M6x10 BLACK     | 2994127       | 2          |
| 3B       | X005544         | 1912111       | 2          |
| 3C       | M16 type S      | 1911653       | 2          |
| 3D       | 125cm           | 1004094       | 2          |
| 3E       | D22             | 1994012       | 2          |
| 3F*      | ---             | 1993218       | (x2) 0,95m |
| 3G       | D12x107         | 1905243       | 2          |
| 4        | Wheel 24" (PU)  | 1916956 (PU)  | 1          |
|          | Wheel 24" (LU)  | 1916389 (AIR) | 1          |
| 5**      | Y018538         | 17710         | 2          |
| 6        | D5x18           | 1001589       | 2          |
| 7        | X005513         | 1992600       | 2          |
| 8        | XP02494         | GP17999       | 2          |
| 9        | M5X40 ISO 7379  | 2994002       | 4          |
| 10       | M5 DIN125       | 2994069       | 4          |
| 11       | M5 black        | 2994006       | 4          |
| 12       | M6 washer black | 1001564       | 2          |
| 13       | ---             | 1910149(R)    | 1          |
|          | ---             | 1910148(L)    | 1          |

\* Adjust wire length to the backrest length if needed /  
Dopasowac dlugosc pancerza do wysokosci oparcia

\*\* Use optionally if handle is too short / Uzyc w wypadku  
zbyt krotkiej raczki

EP00392 / B

|                                     |            |            |            |   |
|-------------------------------------|------------|------------|------------|---|
| TRIGO B74/B80 DETACHABLE DRUM BRAKE |            |            |            | Plan No.  |
| WHEELS                              |            |            |            | EP00392/ B  |
| Scale: 1:4                          | Drafted    | Verified   | Confirmed  |  |
| Date                                | 17/04/2019 | 14/07/2020 | 14/07/2020 |   |
| Initials                            | FD         | PM         | TK         |   |